

Aurolab Procoat Plus

Precision Viscoelastic Solution

Sodium hyaluronate 3% w/v with Chondroitin Sulphate 4% w/v Ophthalmic Solution



AUROLAB PROCOAT PLUS: ENHANCED PROTECTION FOR ADVANCED SURGICAL PRECISION



- High purity blend for superior support
- Sterile and isotonic formulation
- Viscoelastic properties for controlled maneuvers
- Maintains anterior chamber depth
- Biocompatible and easy to remove
- Wide range of applications
- Preservative-free for maximum safety
- Ready-to-use convenience
- Optimal rheological properties for clear visibility
- Storage and supply

AUROLAB PROCOAT: ENGINEERED FOR EXCEPTIONAL **SURGICAL MANEUVERS**

PROPERTIES

CHOOSE AUROLAB PROCOAT PLUS AND PROCOAT

PROPERTIES Concentration	Sofum hydrorate 38 w/w with Chordroids subrate 48 w/ cytus as a second control of the control of
Sodium Hyaluronate	30mg/mL (3% w/v)
Chondroitin Sulphate	40mg/mL (4% w/v)
Molecular Weight	
Sodium Hyaluronate	>1M Daltons
Chondroitin Sulphate	>20000 Daltons
рН	6.80 to 7.60
Osmolality	200 to 400 mOsm/Kg
Viscosity at 25°C	20,000 to 60,000 cPs
Refractive Index	1.330 to 1.350 at 25°C
Specific Gravity	1.000 to 1.040 at 25°C
Spectral Transmittance	>90%

Concentration	has been been produced as the frequency
Sodium Hyaluronate	20mg/mL (2% w/v)
Chondroitin Sulphate	20mg/mL (2% w/v)
Molecular Weight	
Sodium Hyaluronate	>3M Daltons
Chondroitin Sulphate	>20000 Daltons
pH	6.80 to 7.60
Osmolality	200 to 400 mOsm/Kg
Viscosity at 25°C	20,000 to 60,000 cPs
Refractive Index	1.330 to 1.350 at 25°C
Specific Gravity	1.000 to 1.040 at 25°C
Spectral Transmittance	>90%

Information published in this catalogue is subject to change without notification

For a demo or more information on how your patients can benefit from Aurolab products, call at 1800 103 7321 or email us at info@aurolab.com

AUPOLAB